

Single Molecule Microscope - TIRF

<https://search.researchequipment.wur.nl/SearchDetail.aspx?deviceid=11ac8726-d59f-4530-91b8-0c79af8d1239>

Brand

home-built

Type



Contact

Johannes Hohlbein (johannes.hohlbein@wur.nl)

Arjen Bader (arjen.bader@wur.nl)

Organisation

Agrotechnology & Food Sciences Group

Department

MICROSPECTR. CENTR. (MSC)

Description

Our single-molecule microscopy (SMM) setup is a home-built total internal reflection fluorescence (TIRF) microscope with 4 (alternating) laser excitation channels and three emission channels. It can therefore be used for single-molecule Förster Resonance Energy Transfer (FRET) spectroscopy, as well as super-resolution microscopy techniques such as photoactivated localization microscopy (PALM) and stochastic optical reconstruction microscopy (STORM).

Principle of TIRF

In TIRF microscopy the laser light is focused into the rim of the backfocal plane of the microscope objective such as that the laser light is totally reflected at the glass/ water interface. An evanescent field is created in which intensity decays within a few hundred nanometre. Thereby, only a small volume above the coverslide is illuminated.

smFRET is a well-established technique to measure distances between two fluorophores with high spatial (nanometre) and temporal (millisecond) resolution. In FRET, the proportion of transferred energy between donor fluorophore and an acceptor fluorophore depends strongly on the distance r between them. This makes the technique useful to study, for example, conformational dynamics of enzymes such as DNA polymerases. In the SMM set-up, we detect fluorophores close to the surface within an area of roughly 50 by 50 micrometer by imaging the emitted with an emCCD camera. We can detect hundreds of single molecules in parallel with TIRF microscopy.

A very important addition to the concept of single-molecule FRET was the development of ALEX (Alternating Laser Excitation): instead of exciting only the donor fluorophore and measuring the energy transfer to the acceptor at all times, the acceptor is repeatedly directly excited in order to probe its existence. If there is, for example, no red fluorescence detectable after red excitation, the acceptor is already bleached and is not participating in FRET any more. The detection of red fluorescence after red excitation allows the introduction of a new parameter called stoichiometry, which relates the green and red fluorescence after green excitation to the overall fluorescence emission after green and red excitation. The figure shows an ES histogram allowing the clear separation of species which cannot easily be resolved on the E histogram alone.

Super-resolution microscopy (SRM)

The SMM system can be used for super-resolution imaging with techniques like stochastic optical reconstruction microscopy (STORM) or photoactivated localization microscopy (PALM). The location of individual single molecules can be determined with high accuracy (typically 20 nm), but when two molecules are less than ~200 nm apart they cannot be distinguished (i.e. the resolution is not sufficient). Super-resolution images can be recorded by ensuring that most of the fluorescent labels are in an off-state for most of the time. A sequence of images is now recorded where in each image only a few isolated single molecules emit fluorescence stochastically (STORM) or after photoactivation (PALM). This way the locations of the individual molecules are determined and an image of all these locations can be composed.

Technical Details

- TIRF microscope with 4 excitation wavelengths (405, 473, 561 and 642 nm). Alternating laser excitation with millisecond time resolution.
- Detection by emCCD (Andor iXon) with three channels
- for GFP / fluorescein / Alexa488 / Atto488 (473 nm excitation)
- for DsRED / mCherry / rhodamine / Alexa568 / Atto565 (561 nm excitation)
- for Cy5 / Alexa647 / Atto647 (642 nm excitation)
- Up to 60 frames per second
- Programmable alternating laser excitation schemes (including PALM)

Applications

- Single molecule FRET is used to measure molecular scale distances (2-10 nm) between fluorescence tags in individual proteins / protein complexes / protein-DNA complexes. Also, it can follow binding dynamics.
- Super-resolution microscopy can resolve structures in cells with a resolution of 20 nm. The proteins of interest are selectively imaged by special fluorescence probes.

Complementary Techniques

Super-resolution microscope for PALM and STORM

Publications

smTIRF PPT presentation, Johannes Hohlbein, , http://www.wageningenur.nl/upload_mm/1/d/0/8f66a338-6765-4ac2-abd0-e22dafbda8be_smTIRF%20-%20Johannes%20Hohlbein%2013062013.pdf

Camera-based single-molecule FRET detection with improved time resolution, S. Farooq and J. Hohlbein, Physical Chemistry Chemical Physics, <http://pubs.rsc.org/en/content/articlelanding/2015/cp/c5cp04137f#!divAbstract>

Blog Johannes Hohlbein, , , <http://www.jhohlbein.com/>